

# Microscopy and Staining

## Introduction

Development of Microscopy: Anton Van Leeuwenhoek (1632-1723)

“In the year of 1657 I discovered very small living creatures in rain water”


## Principles of Microscopy

### Metric Units:

Micrometer ( $\mu\text{m}$ )  $1\mu\text{m} = 10^{-6}\text{ m}$  (or 0.000039 in)

Nanometer (nm)  $1\text{nm} = 10^{-9}\text{ m}$  (or 0.000000039 in)

Angstrom (Å)  $1\text{Å} = 10^{-10}\text{ m}$  (or 0.0000000039 in)



Most prokaryotes (Bacteria) range from 0.5 – 2.0  $\mu\text{m}$  (1-10  $\mu\text{m}$ )

Viruses range from 10 nm – 100nm


Yeast cells: 10  $\mu\text{m}$

## Light Microscopy

- The compound light microscope
- Parfocal Microscope
- Dark-Field Microscopy
- Phase-Contrast Microscopy
- Electron Microscopy



# Anatomy of a Microscope



© OnFocus Laboratories

## Techniques of Light Microscopy

- Wet Mounts





- Smears

- Heat fixing

- Simple Staining

- Differential Staining: Gram Stain

- Special Staining: i.e. Endospore staining




## Differential Staining – Gram Staining

Developed by Hans Christian Gram in 1884

Gram Positive – Bacteria stain purple by crystal violet

Gram Negative – Bacteria stain red by safranin and do not retain crystal violet stain



## Protocol

### Materials Needed:

- 1.) Clean glass microscope slide
- 2.) Sterile inoculating loop
- 3.) Water dropper or pipette
- 4.) Beaker of water
- 5.) Bunsen burner
- 6.) Forceps
- 7.) Marker or wax pencil
- 8.) Culture of bacteria
- 9.) Gram stain solutions


### Slide Preparation:

- Obtain a clean glass slide.
- Draw a circle on the center of the slide where you will place your bacterial culture.
- Place a droplet of water into the circle.
- Remove one colony of bacterial culture using a sterile inoculating loop and smear the bacteria into the water droplet.
- Holding the end of the slide with forceps, heat fix the bacteria by passing the bottom of the slide through the Bunsen burner flame for 5 seconds.
- Let the slide cool before staining.

### Gram Staining Procedure:

- Flood the smear with crystal violet (primary stain) until is completely covered. Leave it on for 30 seconds and then discard into the sink.
- Rinse the slide with water from a wash bottle of distilled water. Drain off excess water.

- Flood the slide with Gram's iodine solution (mordant) and let it stand for 30 seconds and wash with water.
- Hold slide at 45° angle and slowly add 95% ethyl alcohol over for 15 seconds, until no more crystal violet runs off, (timing is critical in this step). Immediately wash with water for few seconds.
- Cover the smear with safranin (counter stain) for 30 seconds, and wash with water.
- Blot the slide carefully with bibulous paper, and air dry.
- Examine the slide using the phase contrast microscope using immersion oil.
- Describe the color and morphology of your stained bacteria:



---

---

---

---

## Yeast Study

If time permits,

*Saccharomyces cerevisiae* will be examined.

- Living culture: Place two or three loopfuls of *S. cerevisiae* on the slide with a drop of methylene blue stain. Place a cover slip and examine under high-dry or immersion oil.
- Look for the shape of the cells, nucleus, vacuoles, blastospore or bud.